

SwabExpress: An end-to-end protocol for extraction-free COVID-19 testing

Sanjay Srivatsan^{*1}, Sarah Heidl^{*2}, Brian Pfau^{*2}, Beth K. Martin¹, Peter D. Han², Weizhi Zhong², Katrina van Raay², Evan McDermot², Jordan Opsahl², Luis Gamboa², Nahum Smith², Melissa Truong², Shari Cho², Kaitlyn A. Barrow³, Lucille M. Rich³, Jeremy Stone², Caitlin R. Wolf⁴, Denise J. McCulloch⁴, Ashley E. Kim⁴, Elisabeth Brandstetter², Sarah L. Sohlberg⁴, Misja Ilcisin⁵, Rachel E. Geyer⁶, Wei Chen¹, Jase Gehring¹, Seattle Flu Study Investigators, Sriram Kosuri^{7,8}, Trevor Bedford^{1,2,5}, Mark J. Rieder², Deborah A. Nickerson^{1,2}, Helen Y. Chu^{2,4}, Eric Q. Konnick^{2,9}, Jason S. Debley³, Jay Shendure^{1,2,10}, Christina M. Lockwood^{1,2,9,11} and Lea M. Starita^{1,2,11}

* denotes equal contributors

1. Department of Genome Sciences, University of Washington, Seattle WA, USA

2. Brotman Baty Institute For Precision Medicine, Seattle WA, USA

3. Center for Immunity and Immunotherapies, Seattle Children's Research Institute, Seattle WA, USA

4. Department of Allergy and Infectious Disease, University of Washington, Seattle WA, USA

5. Vaccine and Infectious Disease Division, Fred Hutchinson Cancer Research Center, Seattle, WA, USA

6. Department of Family Medicine, University of Washington, Seattle, Washington, USA.

7. Octant, Inc. Emeryville CA, USA

8. Department of Chemistry and Biochemistry, University of California, Los Angeles, Los Angeles CA, USA

9. Department of Laboratory Medicine and Pathology, Seattle WA, USA

10. Howard Hughes Medical Institute. Seattle WA, USA

11. Correspondence should be addressed to lstarita@uw.edu or tinalock@uw.edu

Lea Starita

University of Washington

Box 355065

Seattle, WA 98195

lstarita@uw.edu

Tel: 206-685-5495

Fax: 206-685-7301

Christina Lockwood

University of Washington

Box 357110

Seattle, WA 98195

tinalock@uw.edu

Tel: 206-598-2142

Fax: 206-543-3644

1

2

344Structured Abstract

4

545

6

746Background

8

947The urgent need for massively scaled clinical testing for SARS-CoV-2, along with global shortages

1048of critical reagents and supplies, has necessitated development of streamlined laboratory testing

1149protocols. Conventional nucleic acid testing for SARS-CoV-2 involves collection of a clinical

1250specimen with a nasopharyngeal swab in transport medium, nucleic acid extraction, and

1351quantitative reverse transcription PCR (RT-qPCR) (1). As testing has scaled across the world, the

1452global supply chain has buckled, rendering testing reagents and materials scarce (2). To address

1553shortages, we developed SwabExpress, an end-to-end protocol developed to employ mass

1654produced anterior nares swabs and bypass the requirement for transport media and nucleic acid

1755extraction.

1856

1957Methods

20

2158We evaluated anterior nares swabs, transported dry and eluted in low-TE buffer as a direct-to-

2259RT-qPCR alternative to extraction-dependent viral transport media. We validated our protocol of

2360using heat treatment for viral inactivation and added a proteinase K digestion step to reduce

2461amplification interference. We tested this protocol across archived and prospectively collected

2562swab specimens to fine-tune test performance.

2663

2764Results

28

2965After optimization, SwabExpress has a low limit of detection at 2-4 molecules/uL, 100% sensitivity,

3066and 99.4% specificity when compared side-by-side with a traditional RT-qPCR protocol employing

3167extraction. On real-world specimens, SwabExpress outperforms an automated extraction system

3268while simultaneously reducing cost and hands-on time.

3369

70 **Conclusion**

71 SwabExpress is a simplified workflow that facilitates scaled testing for COVID-19 without
72 sacrificing test performance. It may serve as a template for the simplification of PCR-based clinical
73 laboratory tests, particularly in times of critical shortages during pandemics.

74

75

76

77

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

78 Introduction

79 Since the first reported cases in the winter of 2019, the spread of the novel beta-coronavirus
80 SARS-CoV-2 has grown into a global pandemic. The virus spreads easily from person to person
81 and is often carried by asymptomatic individuals (3). These viral properties, in conjunction with a
82 lack of an effective centralized response or societal adherence to public health recommendations,
83 has led to a continued persistence of the pandemic throughout the United States (4). It is widely
84 recognized that increased testing capacity can ameliorate the outbreak (5,6), but the prohibitive
85 cost of testing materials and reagents as well as global supply chain problems continue to thwart
86 efforts to reach the required scale.

87

88 Since the beginning of the pandemic, the gold standard for SARS-CoV-2 detection has been RNA
89 extraction followed by reverse-transcription quantitative polymerase chain reaction (RT-qPCR).
90 Specimens are traditionally collected as nasopharyngeal (NP) specimens (3) by healthcare
91 professionals and transported in viral media (e.g. Universal Transport Media (UTM)). Worldwide
92 reliance on this template protocol has led to global shortages in swabs, viral media, and laboratory
93 reagents. These shortages continue to plague testing labs and impede efforts to scale. Prior
94 literature (7,8) and the work of United Health/Quantigen (9), have established that swabs collected
95 without transport media are acceptable for nucleic acid detection-based diagnostics, eliminating
96 the reliance on UTM. Extraction-free protocols have also been developed to remove the need for
97 RNA extraction reagents and streamline testing protocols. Saliva specimens have been shown to
98 be particularly amenable to extraction-free testing protocols. For example, SalivaDirect™ – a
99 protocol for performing SARS-CoV-2 RT-qPCR on saliva specimens without extraction (10) – had
100 a sensitivity of 89% compared to traditionally processed anterior nares (AN) or oropharyngeal
101 (OP) swabs, demonstrating the viability of extraction-free protocols. Unlike saliva, extraction-free
102 methods for nasal swabs have been less sensitive than conventional protocols—likely due to PCR
103 inhibition from transport media or saline (10–16).

104
105 Here we describe the development of an UTM and extraction-free protocol for anterior nasal dry-
106 swabs that is compatible with RT-qPCR and does not sacrifice test performance. This protocol,
107 which we have coined 'SwabExpress', has a low limit of detection, high sensitivity, high specificity,
108 and superior test performance when compared to conventional extraction-based RT-qPCR
109 protocols. We further identify and ameliorate two distinct failure modes for extraction free RT-
110 qPCR-based testing. Widespread adoption of this approach and others like it could result in a
111 dramatic increase in testing capacity, decrease consumables used during testing, and ultimately
112 help curb the spread of SARS-CoV-2.

113

114 **Methods**

115 *Collection of Nasal Swabs* - For preliminary studies, individuals who tested positive for SARS-
116 CoV-2 through clinical testing were identified and recruited into a study of home-based, self-
117 collected home swabs (28). After providing consent, enrolled participants were supplied a Swab
118 and Send kit (19) containing two swabs (Copan FloqSwab 56380CS01) delivered to their home
119 via 2-hour delivery and were provided instructions to self collect two mid-turbinate swabs.
120 Participants placed one swab in a tube with UTM (Becton Dickinson PN 220220) and the other in
121 an empty, dry 15 mL conical tube for transport. For all other studies, anterior nares (US Cotton
122 #3, distributed by Steripack) swabs were collected by the Seattle Flu Study, Husky Coronavirus
123 Testing Program (HCT) (26) or the Seattle Coronavirus Assessment Network (SCAN) (29).
124 Anterior nares swabs were transported in a sterile, empty conical tube directly to the lab by HCT
125 technicians. SCAN swabs were packaged by the participant according to kit instructions and sent
126 to the Brotman Baty Institute / Northwest Genomics Center, utilizing standard International Air
127 Transport Association shipping procedures by courier at ambient temperature. These IRB-
128 supervised studies were public health surveillance programs and enrolled both symptomatic and
129 asymptomatic participants. Informed consent was obtained from adult participants and

1
2
3 130 parents/permanent legal guardians of participant children. Archived and fresh convenience
4
5 131 specimens from these studies were chosen at random for use in the current study.
6
7 132
8
9 133 *Usability Study* - To recruit a sufficient number of children for the prospective usability study,
10
11 134 participants were recruited that met broad eligibility criteria: 1. No COVID-19 symptoms, 2. no
12
13 135 prior self-swab experience and 3. no prior medical or laboratory training. We obtained informed
14
15 136 consent from adult participants and parents/permanent legal guardians of participant children.
16
17 137
18
19 138 *Swab rehydration and elution* - All work was performed within a class II biosafety cabinet with
20
21 139 appropriate precautions. For preliminary studies each mid turbinate dry swab was placed into a
22
23 140 1.5 mL microfuge tube, then cut using a sterile razor blade. Next, 200 μ L of low-TE [10mM Tris-
24
25 141 HCl pH 7.5 (T2319-1L, Sigma), 0.1mM EDTA (15575020, Invitrogen)] was added to each tube
26
27 142 and vortexed for 30 seconds. To test various buffers, 45 μ L of this solution was removed and
28
29 143 added to either 5 μ L of low-TE or 5 μ L of 10% Triton-X (X100-500ML, Sigma Aldrich). These two
30
31 144 specimens constitute the undiluted eluate from the dry swabs.
32
33 145
34
35 146 For all other studies anterior nares swabs were rehydrated in 1 mL low-TE prepared in UltraPure
36
37 147 Water (Life Technologies PN 10977023). Specimens were vortexed for 30 seconds or shaken for
38
39 148 1 minute and allowed to incubate at room temperature for at least 10 minutes before transfer to
40
41 149 Matrix tubes (Thermo Fisher).
42
43 150
44
45 151 *RNA extraction of specimens* - 200 μ L of eluate was extracted on the Magna Pure 96 using a
46
47 152 DNA and Viral NA Small Volume Kit (Roche, 06543588001) with the universal small volume
48
49 153 protocol and eluted into 50 μ L proprietary elution buffer. Or 200 μ L of eluted anterior nares
50
51 154 specimens were extracted on the KingFisher Flex using the MagMAX Viral Pathogen II Nucleic
52
53 155 Acid Isolation Kit with MagMAX™ Viral/Pathogen Ultra Enzyme Mix (Thermo Fisher A48383 and
54
55
56
57
58
59
60

156 A42366) and eluted in 50 μ L (although circa 35 μ L is eluted).

157

158 *SwabExpress Specimen Preparation* - 50 μ L of 94 specimens were transferred to a LoBind 96
159 well plate (Eppendorf 30129512) using a manual 96-well pipetting system (Rainin Liquidator) with
160 low retention tips (Rainin 17014402) with or without 5 μ L of Proteinase K (Thermo Fisher A42363,
161 proprietary concentration). The plate was sealed with foil (Eppendorf 0030127854 and
162 5392000013). Specimens with Proteinase K were incubated at 37°C for 15 minutes in a
163 convection oven (Across International 0853924003042) and then transferred to a second oven
164 for heat inactivation at 95°C for 15 minutes. Specimens without Proteinase K were heat
165 inactivated at 95°C for 30 minutes.

166 *RT-qPCR* - Each RT-qPCR reaction was performed at a final volume of 10 μ L and containing 1X
167 TaqPath RT-qPCR MasterMix (PN A15300, Life Technologies), 0.125X RNase P TaqMan VIC
168 assay (A30064, Life Technologies) or 1X RNase P HEX assay (IDT), 1X SARS-Cov-2 ORF1b
169 FAM assay (PN 4332079, Life Technologies assay# APGZJKF) or 1X Spike (S) gene (PN
170 4332079, Life Technologies assay# APXGVC4) and nuclease-free water (1907076, Thermo
171 Fisher). 5 μ L of specimen was added to each well. Primer sequences were designed against
172 Wuhan-Hu-1 sequence (MN908947.3) and are proprietary to Thermo Fisher. Plates were sealed
173 using optically clear microseal B (Biorad). Each assay was performed in technical duplicate for a
174 total of four RT-qPCR wells per sample. RT-qPCR was then performed on the Applied Biosystems
175 QuantStudio 6 Pro (25°C for 2 minutes, 50°C for 15 minutes, 98°C for 3 minutes, followed by 40
176 cycles of 98°C for 3 seconds and 60°C for 30 seconds). Reported Ct values were obtained from
177 the onboard analysis using predetermined cycle thresholds. Positive controls contained purified
178 nucleic acid with sequence that was amplified by the ORF1b and Spike gene assays.

179

180 The RT-qPCR reaction for the CDC COVID-19 diagnostic test was performed at a final volume of
181 20 μ L. Reactions contained 1X TaqPath RT-qPCR MasterMix , nCOV-N1 FAM or nCOV-N2 FAM

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

182 primer and probe mix (10006713, IDT) and nuclease-free water (1907076, Thermo Fisher) and 5
183 μ L of specimen was added to each well. RT-qPCR was then performed on the QuantStudio 6
184 Pro as above. Reported Ct values were obtained from the onboard analysis using the auto-
185 determined cycle thresholds. Data were analyzed using Excel and R 3.5.

186

187 *Preparation of inactivated viral controls* - Contrived SARS-CoV-2 positive swabs were generated
188 by collecting clinical matrix from a confirmed healthy volunteer and loaded with 2 μ L of diluted
189 heat-inactivated virion (VR-1986HK [1.6e6 virion/ μ L], ATCC).

190

191 *Viral Inactivation Studies* - Viral inactivation studies were performed at the Seattle Children's
192 Research Institute biosafety level 3 facility. 25 μ L of viral stock (isolate USA-WA1/2020 obtained
193 from ATCC BEI Resources) with a titer of 5.8×10^6 pfu/mL was incubated in 200 μ L of TE or TE +
194 0.25% Triton for 10 minutes at room temperature, or in TE at 65°C for 10 minutes. Untreated and
195 treated SARS-CoV-2 was then added neat and at 10-fold dilutions through 10^{-7} to confluent
196 cultures of Vero E6 cells (CRL-1586, ATCC) and 48 hours later cytopathic effects were scored
197 after staining with crystal violet. RNA was isolated from Vero cells using a TRIzol Plus RNA
198 Purification Kit (ThermoFisher) and the amount of SARS-CoV-2 was quantified by RT-qPCR.

199

200 *Retrospective Comparison Studies* - Remnant participant specimens were stored either at 4°C or
201 -80°C and prepared for RT-qPCR by extraction or heat treatment or SwabExpress digestion as
202 described above. Technicians performing testing and clinical directors interpreting results were
203 both blinded to previous test results.

204

205 *Prospective Comparison Studies* - freshly acquired specimens from the SCAN and HCT studies
206 were prepared by extraction or heat treatment or SwabExpress and tested by RT-qPCR in
207 parallel. For prospective analyses, both technicians and clinical directors performed testing and

interpretation blinded to results from the comparator method.

Results

Usability and Reliability of Anterior Nares (AN) Swabs for At-Home Specimen Collection

We first explored the use of anterior nares (AN) swabs for specimen collection. For mass testing purposes, a swab that is widely available, inexpensive, easy to manufacture, and simple for self-collection is critical. The US Cotton #3 swabs fit these specifications; a polyester AN swab that resembles consumer-brand Q-tips (17). For the purposes of scaled observed or at-home self specimen collection or specimen collection for a child, swabbing the anterior nares anatomical site would be more comfortable, accessible and easier to describe to test users leading to fewer mistakes and better specimen collection (18,19).

Therefore, we conducted a usability study to determine both the accuracy and ease of AN swabs in a Swab-and-Send program where at-home specimen collection kits were delivered to participant residences, the participants swabbed themselves or a child while being virtually monitored by clinical study coordinators and then packaged the specimen for return to the molecular testing lab (19,20). After using the specimen collection kit, study participants completed a survey reporting their level of confidence, the kit's ease of use, and the level of discomfort experienced during swabbing. Participants were recruited from the greater Seattle area and spanned a range of ages, races, household income, and educational attainment (online **Supplemental Fig. S1; Supplemental Table S1A-D**).

The results of the usability study were very encouraging. The majority of participants reported only mild-discomfort during specimen collection with 40% of participants reporting no discomfort at all (**Fig. 1A**). A majority of study participants also found the instructions clear and felt confident that they had correctly collected their specimen (**Fig. 1B**). This was confirmed by low observed

1
2
3 234 rates of error during specimen collection using the AN swabs and during packaging for return
4
5 235 (online **Supplemental Table S2A-B**). Molecular testing performed on these self-collected
6
7 236 specimens confirmed this; RT-qPCR detected the human marker RNase P mRNA in 100% of
8
9 237 swabs with an average crossing threshold value of 23.5 (SD 1.7). The amount of RNase P mRNA
10
11 238 recovered from the AN swabs was higher than for unsupervised collection of mid-turbinate swabs,
12
13 239 which had an average crossing threshold (Ct) value of 26.9 (SD 2.5) (**Fig. 1C**). Together these
14
15 240 data indicate that the use of widely available polyester swabs in the anterior nares is a viable and
16
17 241 preferable alternative for at home specimen collection.
18
19
20 242

21
22 243 Handling Dry Swabs in the Clinical Laboratory
23

24 244 Standard viral media such as Universal Transport Medium (UTM) (e.g. COPAN Diagnostics) have
25
26 245 been in short supply over the course of the pandemic. These salt-rich media inhibit direct RT-
27
28 246 qPCR, making RNA extraction a necessity and thus create an additional bottleneck in the testing
29
30 247 process. Furthermore, automated extraction systems are expensive and their reagents and
31
32 248 consumables are also subject to global shortages. Therefore, we focused on eliminating UTM
33
34 249 and extraction from our testing platform. To bypass UTM, we adopted a dry swab transport and
35
36 250 rehydration method validated by Quantigen that has been explored by other clinical testing
37
38 251 laboratories (15,21). Next, to eliminate RNA extraction and enable direct RT-qPCR, we tested
39
40 252 rehydration solutions for their ability to elute contrived SARS-CoV-2 specimens, compatibility with
41
42 253 direct RT-qPCR, and simplicity. We determined that elution in low-TE (10 mM Tris pH 7.5, 0.1
43
44 254 mM EDTA) without other detergents was best suited for direct RT-qPCR (online **Supplemental**
45
46 255 **Fig. S2**). Unlike UTM and other saline solutions, the low ionic strength of low-TE does not inhibit
47
48 256 PCR amplification. Moreover, low-TE can be quickly prepared using reagents commonly found in
49
50 257 laboratories.
51
52 258

Bypassing nucleic acids extraction poses another problem; instead of the virus being inactivated by the denaturing agents during nucleic acid extraction, the specimen eluted from the swab remains potentially infectious for SARS-CoV-2 or other pathogens and poses a risk to laboratory staff. Accordingly, specimens from both conventional UTM and rehydrated dry swabs are processed inside a class II biosafety (BSL-2) cabinet, in accordance with federal regulatory guidance. However, it is practical and beneficial for downstream steps (like preparing RT-qPCR reactions) to take place on a BSL-2 designated bench. Therefore, we compared several inactivation methods to determine which would be easiest without inhibiting PCR or causing a loss of sensitivity. Viral inactivation of coronaviruses can be achieved through the use of either detergent or heat (22). Our prior results demonstrate the negative impact of detergents on RT-qPCR (online **Supplemental Fig. S2**); therefore, we opted to deploy heat inactivation (online **Supplemental Fig. S3**). We used a protocol to heat inactivate at higher temperatures (95°C) for 30 minutes to increase the safety margins. We also determined that this high-heat protocol had the added benefit of stabilizing the sample over time, a result concordant with another SARS-CoV-2 testing protocol in saliva (23).

Performance of Extraction-free RT-qPCR

Having developed an extraction-free RT-qPCR protocol (EF-RT-qPCR), we set out to determine its performance on both contrived and clinical specimens. To assess analytical sensitivity, we first determined this assay's limit of detection (LoD), the minimum number of SARS-CoV-2 RNA molecules that could be detected in greater than 95% of RT-qPCR reactions. To generate these contrived specimens, we inoculated AN swabs with clinical matrix collected from a healthy volunteer with dilutions of heat inactivated SARS-CoV-2. These experiments determined the EF-RT-qPCR analytical sensitivity to be 2 molecules/μl of eluate for the Orf1b assay and 4 molecules/μl of eluate for the S-gene (Spike gene) assay (online **Supplemental Table S3**). This

1
2
3 284 LoD is comparable to the LoD of many other RT-qPCR-based tests that have been issued
4
5 285 Emergency Use Authorization from the FDA (24).
6
7 286
8
9 287 Next we tested the performance of EF-RT-qPCR compared to our clinically validated RT-qPCR
10
11 288 laboratory-developed test on archived AN specimens. In this assay each sample is tested in 4
12
13 289 independent RT-qPCR reactions, comprising two SARS-CoV-2 assays (Orf1b and Spike) in
14
15 290 duplicate, and is multiplexed with a RNase P assay in every well (**Fig. 2A** and online
16
17 291 **Supplemental Fig S4**). Following RT-qPCR, a clinical result is determined by the number of
18
19 292 replicates displaying SARS-CoV-2 amplification: positive (3 or 4 of 4 wells), low-
20
21 293 positive/inconclusive (2 of 4 wells) and negative (0 or 1 of 4 wells). Head-to-head comparison
22
23 294 between EF-RT-qPCR and a reference standard extraction-based RT-qPCR assay on matched
24
25 295 specimens established that EF-RT-qPCR was 100% specific (56/56 negative specimens) and
26
27 296 91.0% sensitive (61/67 - 56 positive and 5 low-positive) (**Fig. 2B**). Comparison of the mean delta
28
29 297 Ct (Δ Ct) values between the two assays showed that eliminating extraction did decrease
30
31 298 analytical sensitivity. We observed an average increase of 1.96, 2.45, and 4.00 cycles for Orf1b,
32
33 299 Spike and RNase P assays, respectively. Indeed, the 6 specimens not detected by EF-RT-qPCR
34
35 300 had an average Ct with the extraction-based RT-qPCR assay of 34.13 for Orf1b and 35.29 for
36
37 301 Spike.
38
39 302
40
41 303 Owing to an unstable supply chain, while validating the EF-RT-qPCR protocol, our clinical
42
43 304 laboratory was forced to switch from the Roche Magna Pure 96 to the Thermo Fisher KingFisher
44
45 305 Flex automated nucleic acids extraction platform. The relative sensitivity, specificity and Δ Ct
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

values between 619 prospective specimens run in parallel on both the KingFisher Flex (extraction) and EF-RT-qPCR were comparable to results of the retrospective study on stored specimens. EF-RT-qPCR detected SARS-CoV-2 in 100% of specimens that were positive by the extraction method with a 99.4% specificity (**Tables S5-S7**).

Addition of Proteinase K reduces amplification interference

After deploying EF-RT-qPCR as our clinical testing platform, we repeatedly observed two undesirable outcomes that were not observed in our validation studies. First, for 0.9% of specimens (n = 383/43,539), amplification of the human RNase P internal control was undetected in 2 or more of the 4 reactions (**Fig. 2B**; online **Supplemental Table S7**). These specimens were classified as “failures” and each test was repeated before releasing the result. Second, for 0.5% of specimens (229/43,539), we sporadically observed the presence of strong amplification (Ct < 30) in a single well for one of the SARS-CoV-2 targets in specimens where the three other wells were undetected (online **Supplemental Table S8**). However, upon repeat RT-qPCR, both with and without extraction, all 4 wells of the SARS-CoV-2 reactions for these specimens were undetected.

We noted that some of the specimens that produced these problematic outcomes had excess mucous or other nasal secretions. Therefore, we hypothesized that the addition of proteinase K (ProK) digestion could ameliorate both RNase P failures and the spurious SARS-CoV-2 amplification by digesting mucins and other potentially interfering proteins in the nasal specimens (25). We compared RT-qPCR results for 1,222 clinical specimens prepared by the 30 minute 95°C heat treatment with those digested with ProK for 15 minutes prior to heat treatment at 95°C for 15 minutes. We observed approximately 10-fold fewer RT-qPCR reactions with failed RNase P

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

amplification -- 27 of 4,888 without ProK vs 2 of 4,888 with ProK -- reducing the failure rate to 0.04% (**Fig. 3A**, online **Supplemental Table S9**), and improved RNase P detection (ΔCt -0.88) (online **Supplemental Fig. S5**). Furthermore, the addition of a ProK digestion step eliminated spurious amplification of SARS-CoV-2 targets.

In the 4,888 specimens processed both with and without ProK, ProK-treated specimens had decreased Ct values (mean decrease of 1.22 for Orf1B, and 0.97 for Spike). This increased sensitivity was also reflected in the ability to accurately classify archived SARS-CoV-2 positive specimens with Ct values > 28 (**Fig. 3B**). Repeatability and reproducibility were also improved with the addition of ProK (**Fig. 3C**). Upon addition of ProK, on SARS-CoV-2 positive samples, our protocol had a higher concordance (93.3%) versus without ProK (90%) or specimens extracted on the KingFisher Flex (86.6%) (online **Supplemental Table S10**). After this optimization we named our final protocol “SwabExpress” -- consisting of a dry AN swab, followed by ProK digestion and direct RT-PCR. Finally, we prospectively compared performance on 1169 specimens run in parallel on the SwabExpress and KingFisher Flex (extraction) platforms. Positive and negative clinical concordance was excellent; there was 100% concordance for positives results, 99.91% concordance across negatives with a small ΔCt value of 0.37 for the Orf1b target and 1.46 for the S target between the two assays (**Fig. 3D**, online **Supplemental Table S11**).

SwabExpress is compatible with other SARS-CoV-2 RT-qPCR assays

Our laboratory-developed test uses custom Orf1b and Spike-gene assays for detecting SARS-CoV-2. To establish that the SwabExpress protocol was compatible with the widely used CDC N1 and N2 assays, we performed RT-qPCR on 75 positive specimens and 92 negative specimens

with the N1 and N2 assays performed in parallel on the SwabExpress platform and extraction-based RT-qPCR platform. The results were 100% concordant between our custom assays and the CDC assays. Ct values for positive samples were delayed when prepared by SwabExpress protocol compared to the Roche Magna Pure 96. However, this difference did not change the clinical interpretation of these samples (**fig. S6**). For the N1 assay, extracted specimens had an average Ct of 19.22 ± 3.67 versus 21.79 ± 4.33 with SwabExpress (Δ Ct of 2.57). For the N2 assay, extracted specimens had an average Ct of 18.31 ± 3.73 versus Cts of 19.80 ± 3.72 for SwabExpress (Δ Ct of 1.49) (**Table S12**).

SwabExpress is time and cost effective

A dry-swab, extraction-free RT-qPCR protocol comprises the minimal components of a diagnostic test. Although the addition of a proteinase K digestion adds \$0.14 USD to the reagent cost for each sample, this cost is warranted. The addition of proteinase K reduces the repeat rate, reduces the chances of a false positive result from interfering substances during PCR amplification, improves the performance of the test, and in our hands, outperformed a suboptimal yet widely used automated extraction system (Thermo KingFisher Flex™).

Upon adoption, SwabExpress approximately doubled laboratory capacity. First, hands-on technician time previously spent preparing and running extraction systems, went towards accessioning and processing additional samples. Second, the SwabExpress protocol increases scale by using a convection oven that can process up to 6 96-well plates simultaneously. This throughput greatly exceeds the single 96-well plate processed by commercial automated extraction systems. Further scaling of the SwabExpress protocol can be accomplished through

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

the purchase of additional or larger ovens, although RT-qPCR instruments used during amplification and readout still pose a substantial bottleneck in the testing protocol.

Along with the substantial cost of purchasing automated extractors, the consumables required for their operation cost between \$4 and \$5 per sample. By eliminating extraction, and transport medium, SwabExpress reduces the associated costs by more than 90% (~\$0.20 per sample). In all, SwabExpress offers a time and cost-saving alternative to nucleic acids extraction using readily available reagents, which reduces dependence on a heavily burdened supply chain (**Fig. 4, Table S13**).

Discussion

Here we present SwabExpress, an end-to-end diagnostic platform optimized for faster and simpler low-cost detection of SARS-CoV-2 from nasal swabs without the use of nucleic acid extraction (**Figure 4**). This protocol was so named for its ease, rapid turn around and simplicity - dry swabs, without extraction, enhanced with proK digestion. By eliminating transport media and extraction from the workflow, we have decreased cost per sample and reduced supply chain pressure for the lab. Because of the reduced cost and the ability to process many more specimens in parallel, our lab's capacity markedly increased with its adoption. Importantly, we gained efficiency without sacrificing accuracy; our results suggest that the simplified SwabExpress protocol (direct elution from dry swab into low-TE + proteinase K → RT-qPCR) is as sensitive as the conventional PCR protocol (swab → UTM → RNA extraction → RT-qPCR). SwabExpress has supported scaled testing in our lab with over 91,000 tests performed to date and allowed us to

support large testing endeavors such as the Husky Coronavirus Testing Program for the University of Washington (26).

There are some caveats to consider. Even with the addition of proteinase K, specimens with excess mucous fail to amplify RNase P. Since adding this proteinase digestion step, 18/12,991 specimens have had two or more RNase P reactions fail (0.1%) and our laboratory reflexes these few specimens to an extraction protocol. However, 0.1% compares favorably when compared to a protocol with extraction where the failure rate due to failed RNase P is 1% (215/22,546). In addition, the unknown presence of inhibitors precludes comparison of Ct values between specimens; therefore, studies directly comparing Ct values from different specimens may not yield accurate results.

We have observed a marked loss of viral RNA after freeze-thaw cycles for specimens stored in low-TE compared to specimens stored in commercial UTM. We detect a Δ Ct of about 2.5 for specimens after -80°C storage in low-TE, whereas the Δ Ct for specimens retested after storage at -80°C in UTM has historically been negligible. This impacts the ability to use these specimens for downstream applications such as genomic sequencing.

Several improvements can be incorporated into the SwabExpress platform. First, the ability to detect multiple pathogens from one assay can be explored. It is likely that SwabExpress will be compatible with other enveloped viruses such as influenza and respiratory syncytial virus. Multiple targets can be detected in many qPCR systems and the reemergence of these viruses as COVID-19 prevention control measures are relaxed will be of interest to monitor. Second, the most labor intensive part of the SwabExpress protocol is sample accessioning. Receiving individual 10 ml

1
2
3 423 tubes and transferring the eluate to 96-well format takes approximately 2.5 minutes per sample.
4
5 424 Receiving nasal swabs in 96-well compatible, lab-ready transport tubes would streamline the
6
7 425 process considerably (27). Third, incubation times for proteinase K digestion, heat inactivation,
8
9 426 and RT-qPCR could be further optimized to save additional time during the testing process.
10

11 427
12
13 428 Massive scaling and deployment of SARS-CoV-2 testing is essential to curtailing the COVID-19
14
15 429 pandemic, and will likely be necessary well into the future. The protocol evaluated here, including
16
17 430 thousands of real-world, self-collected nasal swabs, would markedly simplify the workflow for RT-
18
19 431 qPCR, the most widely deployed testing paradigm, by eliminating the need for viral transport
20
21 432 media and RNA extraction, both of which are currently experiencing significant supply chain
22
23 433 challenges. Looking forward, we envision that nasal swabs -- self-collected into laboratory ready,
24
25 434 barcoded tubes and transported dry -- could potentially serve as a common input to a range of
26
27 435 SARS-CoV-2 nucleic acid tests for public health surveillance applications. This includes gold-
28
29 436 standard tests like RT-qPCR, but also potentially new modalities like Swab-Seq (11). The
30
31 437 operationalization of the mass distribution and return of such lab-ready collection devices is a
32
33 438 significant effort that should begin now.
34
35

36 439
37
38 440 **Ethics Approval:** Sequencing and analysis of specimens from the Seattle Flu Study, the
39
40 441 Hospitalized and Ambulatory Adults with Respiratory Viral Infections (HAARVI) study and the
41
42 442 SCAN study were approved by the Institutional Review Board at the University of Washington
43
44 443 (protocols STUDY00006181, STUDY00000959, STUDY00007628, STUDY00010432,
45
46 444 STUDY00011148). Informed consent was obtained for all participant specimens.
47
48

49 445
50
51 446 **Acknowledgments:** We would like to thank the Seattle Flu Study, HAARVI, SCAN and Husky
52
53 447 Coronavirus Testing program study participants for their invaluable contributions to this research.
54
55 448 We thank Sydney Floth, Jefferson Nguyen, Ashley Gate, and Gift Nwanne for help with the
56
57
58
59
60

usability study, Catherine Moore of Public Health Wales and Dan Wattendorf, Emily Turner and Karen Heichman of the Bill and Melinda Gates Foundation for helpful advice.

Author Contributions: *All authors confirmed they have contributed to the intellectual content of this paper and have met the following 4 requirements: (a) significant contributions to the conception and design, acquisition of data, or analysis and interpretation of data; (b) drafting or revising the article for intellectual content; (c) final approval of the published article; and (d) agreement to be accountable for all aspects of the article thus ensuring that questions related to the accuracy or integrity of any part of the article are appropriately investigated and resolved.*

Authors' Disclosures or Potential Conflicts of Interest: *Upon manuscript submission, all authors completed the author disclosure form. Disclosures and/or potential conflicts of interest:*

Employment or Leadership: S. Kosuri is an employee of Octant Inc., who have developed SwabSeq, a method that will likely benefit from dry swab and extraction free protocols becoming standard. E.Q. Konnick, Association for Molecular Pathology, College of American Pathologists, Washington State Society of Pathologists.

Consultant or Advisory Role: H.Y. Chu reported consulting with Ellume, Pfizer, The Bill and Melinda Gates Foundation, Glaxo Smith Kline, and Merck. J. Shendure is a consultant with Guardant Health, Maze Therapeutics, Camp4 Therapeutics, Nanostring, Phase Genomics, Adaptive Biotechnologies, and Stratos Genomics.

Stock Ownership: S. Kosuri, Octant Inc.

Honoraria: L.M. Starita, AMP.

Research Funding: The Seattle Flu Study and SCAN are administered by the Brotman Baty Institute for Precision Medicine and funded by Gates Ventures, the private office of Bill Gates. The funder was not involved in the design of the study and does not have any ownership over the management and conduct of the study, the data, or the rights to publish. Funding for the HAARVI study comes from DARPA HR001117S0019, the Bill and Melinda Gates Foundation, Emergent Ventures grant to H.Y. Chu. L.M. Starita, and J. Shendure are funded by 1RM1HG010461-01 from the NHGRI and J. Shendure is an Investigator of the Howard Hughes Medical Institute. J.S. Debley is funded by NIH/NIAID K24AI150991-01S1, which supported the viral inactivation studies. Additional funding for this project came from the University of Washington in support of the Husky Coronavirus Testing Program with funds from the United States Senate and House of Representative Bill 748, Coronavirus Aid, Relief, and Economic Security Act (CARES Act). H.Y. Chu has received research funding from Gates Ventures, Sanofi Pasteur, and support and reagents from Ellume and Cepheid outside of the submitted work. J. Shendure has a research collaboration with Illumina. S. Kosuri, Testing For America supporting separate work on SwabSeq development, Ginkgo Bioworks supporting separate work on CLIA lab formation.

Expert Testimony: S. Kosuri, Broad Institute.

Patents: S. Kosuri, Ginkgo Bioworks - Licensing for SwabSeq multiplexed assay development, Harvard U - Royalties from patents relating to genome editing, UCLA - Royalties from patents relating to multiplexed assays.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

Role of Sponsor: The funding organizations played no role in the design of study, choice of enrolled patients, review and interpretation of data, preparation of manuscript, or final approval of manuscript.

Seattle Flu Study Investigators Principal Investigators: Helen Y. Chu^{1,7}, Michael Boeckh^{1,2,7}, Janet A. Englund^{3,7}, Michael Famulare⁴, Christina M. Lockwood^{1,7}, Barry R. Lutz^{5,7}, Deborah A. Nickerson^{6,7}, Mark J. Rieder⁷, Lea M. Starita^{6,7}, Matthew Thompson⁹, Cécile Viboud¹⁰, Jay Shendure^{6,7,8}, and Trevor Bedford^{2,6,7}

Affiliations:

- 1 Department of Medicine, University of Washington
- 2 Vaccine and Infectious Disease Division, Fred Hutchinson Cancer Research Center
- 3 Seattle Children’s Research Institute
- 4 Institute for Disease Modeling
- 5 Department of Bioengineering, University of Washington
- 6 Department of Genome Sciences, University of Washington
- 7 Brotman Baty Institute for Precision Medicine
- 8 Howard Hughes Medical Institute
- 9 Department of Family Medicine, University of Washington
- 10 National Institutes of Health

References

1. Cheng MP, Papenburg J, Desjardins M, Kanjilal S, Quach C, Libman M, et al. Diagnostic Testing for Severe Acute Respiratory Syndrome–Related Coronavirus 2 [Internet]. *Annals of Internal Medicine*. 2020. page 726–34. Available from: <http://dx.doi.org/10.7326/m20-1301>
2. Botti-Lodovico Y, Rosenberg E, Sabeti PC. Testing in a Pandemic — Improving Access, Coordination, and Prioritization [Internet]. *New England Journal of Medicine*. 2021. page 197–9. Available from: <http://dx.doi.org/10.1056/nejmp2025173>
3. Cevik M, Tate M, Lloyd O, Maraolo AE, Schafers J, Ho A. SARS-CoV-2, SARS-CoV, and MERS-CoV viral load dynamics, duration of viral shedding, and infectiousness: a systematic review and meta-analysis. *Lancet Microbe*. Elsevier BV; 2021;2:e13–22.
4. Brauner JM, Mindermann S, Sharma M, Johnston D, Salvatier J, Gavenčiak T, et al. Inferring the effectiveness of government interventions against COVID-19. *Science* [Internet]. American Association for the Advancement of Science; 2021 [cited 2021 Apr 23];371. Available from: <https://science.sciencemag.org/content/371/6531/eabd9338.abstract>
5. Larremore DB, Wilder B, Lester E, Shehata S, Burke JM, Hay JA, et al. Test sensitivity is secondary to frequency and turnaround time for COVID-19 surveillance. *medRxiv* [Internet]. 2020; Available from: <http://dx.doi.org/10.1101/2020.06.22.20136309>
6. Hellewell J, Russell TW, Beale R, Kelly G, Houlihan C, Nastouli E, et al. Estimating the effectiveness of routine asymptomatic PCR testing at different frequencies for the detection of SARS-CoV-2 infections [Internet]. *bioRxiv*. medRxiv; 2020. Available from: <http://medrxiv.org/lookup/doi/10.1101/2020.11.24.20229948>
7. Moore C, Corden S, Sinha J, Jones R. Dry cotton or flocked respiratory swabs as a simple collection technique for the molecular detection of respiratory viruses using real-time NASBA [Internet]. *Journal of Virological Methods*. 2008. page 84–9. Available from: <http://dx.doi.org/10.1016/j.jviromet.2008.08.001>
8. Emerson J, Cochrane E, McNamara S, Kuypers J, Gibson RL, Campbell AP. Home Self-Collection of Nasal Swabs for Diagnosis of Acute Respiratory Virus Infections in Children With Cystic Fibrosis. *J Pediatric Infect Dis Soc*. 2013;2:345–51.
9. U.S. Department of Health and Human Services Food and Drug Administration. Policy for Coronavirus Disease-2019 Tests During the Public Health Emergency [Internet]. US Food and Drug Administration. 2020 [cited 2021 Mar 20]. Available from: <https://www.fda.gov/media/135659/download>
10. Vogels CBF, Brackney D, Wang J, Kalinich CC, Ott I, Kudo E, et al. SalivaDirect: Simple and sensitive molecular diagnostic test for SARS-CoV-2 surveillance [Internet]. *Infectious Diseases (except HIV/AIDS)*. medRxiv; 2020. Available from: <https://www.medrxiv.org/content/10.1101/2020.08.03.20167791v1.abstract>
11. Smyrlaki I, Ekman M, Lentini A, Rufino de Sousa N, Papanicolaou N, Vondracek M, et al. Massive and rapid COVID-19 testing is feasible by extraction-free SARS-CoV-2 RT-PCR [Internet]. *Infectious Diseases (except HIV/AIDS)*. medRxiv; 2020. Available from:

1
2
3 555 <https://www.medrxiv.org/content/10.1101/2020.04.17.20067348v4.abstract>
4
5 556 12. Lübke N, Senff T, Scherger S, Hauka S, Andrée M, Adams O, et al. Extraction-free SARS-
6 557 CoV-2 detection by rapid RT-qPCR universal for all primary respiratory materials. J Clin
7 558 Virol. 2020;130:104579.
8
9 559 13. Bruce EA, Huang M-L, Perchetti GA, Tighe S, Laaguiby P, Hoffman JJ, et al. Direct RT-
10 560 qPCR detection of SARS-CoV-2 RNA from patient nasopharyngeal swabs without an RNA
11 561 extraction step. PLoS Biol. 2020;18:e3000896.
12
13 562 14. Merindol N, Pépin G, Marchand C, Rheault M, Peterson C, Poirier A, et al. SARS-CoV-2
14 563 detection by direct rRT-PCR without RNA extraction. J Clin Virol. 2020;128:104423.
15
16 564 15. Padgett LR, Kennington LA, Ahls CL, Samarasinghe DK, Tu Y-P, Wallander ML, et al.
17 565 Polyester Nasal Swabs Collected in a Dry Tube are a Robust and Inexpensive, Minimal
18 566 Self-Collection Kit for SARS-CoV-2 Testing [Internet]. Available from:
19 567 <http://dx.doi.org/10.1101/2020.10.09.20210302>
20
21 568 16. Ñique AM, Coronado-Marquina F, Mendez Rico JA, García Mendoza MP, Rojas-Serrano
22 569 N, Simas PVM, et al. A faster and less costly alternative for RNA extraction of SARS-CoV-2
23 570 using proteinase k treatment followed by thermal shock. PLoS One. 2021;16:e0248885.
24
25 571 17. U.S. Department of Health and Human Services Food and Drug Administration. FDA,
26 572 Gates Foundation, UnitedHealth Group, Quantigen, and U.S. Cotton Collaborate to
27 573 Address Testing Supply Needs [Internet]. US Food and Drug Administration. 2020 [cited
28 574 2021 Apr 23]. Available from: [https://www.fda.gov/news-events/press-](https://www.fda.gov/news-events/press-announcements/coronavirus-covid-19-update-fda-gates-foundation-unitedhealth-group-quantigen-and-us-cotton)
29 575 [announcements/coronavirus-covid-19-update-fda-gates-foundation-unitedhealth-group-](https://www.fda.gov/news-events/press-announcements/coronavirus-covid-19-update-fda-gates-foundation-unitedhealth-group-quantigen-and-us-cotton)
30 576 [quantigen-and-us-cotton](https://www.fda.gov/news-events/press-announcements/coronavirus-covid-19-update-fda-gates-foundation-unitedhealth-group-quantigen-and-us-cotton)
31
32 577 18. Truong M, Pfau B, McDermot E, Han PD, Brandstetter E, Richardson M, et al. Comparable
33 578 specimen collection from both ends of at-home mid-turbinate swabs. J Clin Microbiol
34 579 [Internet]. 2021; Available from: <http://dx.doi.org/10.1128/JCM.03073-20>
35
36 580 19. Kim AE, Brandstetter E, Wilcox N, Heimonen J, Graham C, Han PD, et al. Evaluating
37 581 Specimen Quality and Results from a Community-Wide, Home-Based Respiratory
38 582 Surveillance Study. J Clin Microbiol [Internet]. 2021; Available from:
39 583 <http://dx.doi.org/10.1128/JCM.02934-20>
40
41 584 20. Chu HY, Englund JA, Starita LM, Famulare M, Brandstetter E, Nickerson DA, et al. Early
42 585 Detection of Covid-19 through a Citywide Pandemic Surveillance Platform. N Engl J Med
43 586 [Internet]. 2020; Available from: <http://dx.doi.org/10.1056/NEJMc2008646>
44
45 587 21. Parikh BA, Wallace MA, McCune BT, Burnham C-AD, Anderson NW. The Effects of “Dry
46 588 Swab” Incubation on SARS-CoV-2 Molecular Testing. [Epub] J Appl Lab Med March 29,
47 589 2021, as doi: 10.1093/jalm/jfab010.
48
49 590 22. Darnell MER, Taylor DR. Evaluation of inactivation methods for severe acute respiratory
50 591 syndrome coronavirus in noncellular blood products. Transfusion 2006;46:1770–7.
51
52 592 23. Ranoa DRE, Holland RL, Alnaji FG, Green KJ, Wang L, Brooke CB, et al. Saliva-Based
53 593 Molecular Testing for SARS-CoV-2 that Bypasses RNA Extraction [Internet]. Cold Spring
54 594 Harbor Laboratory. 2020 [cited 2020 Dec 31]. page 2020.06.18.159434. Available from:

- 595 <https://www.biorxiv.org/content/10.1101/2020.06.18.159434v1>
- 596 24. Resilience Health [Internet]. [cited 2021 Mar 19]. Available from:
597 <https://www.resiliencehealth.com/tests.html>
- 598 25. Mallmann L, Schallenberger K, Demolliner M, Eisen AKA, Hermann BS, Heldt FH, et al.
599 Pre-treatment of the clinical sample with Proteinase K allows detection of SARS-CoV-2 in
600 the absence of RNA extraction [Internet]. Available from:
601 <http://dx.doi.org/10.1101/2020.05.07.083139>
- 602 26. Weil AA, Sohlberg SL, O'Hanlon JA, Casto AM, Emanuels AW, Lo NK, et al. SARS-CoV-2
603 epidemiology on a public university campus in Washington state [Internet]. bioRxiv.
604 medRxiv; 2021. Available from:
605 <http://medrxiv.org/lookup/doi/10.1101/2021.03.15.21253227>
- 606 27. Rhinostics product information [Internet]. Rhinostics. [cited 2021 Jun 21]. Available from:
607 <https://rhinostics.com>
- 608 28. McCulloch DJ, Kim AE, Wilcox NC, Logue JK, Greninger AL, Englund JA, et al.
609 Comparison of Unsupervised Home Self-collected Midnasal Swabs With Clinician-Collected
610 Nasopharyngeal Swabs for Detection of SARS-CoV-2 Infection. JAMA Netw Open.
611 2020;3:e2016382.
- 612 29. Greater Seattle Coronavirus Assessment Network (SCAN) [Internet]. [cited 2021 Feb 19].
613 Available from: <https://scanpublichealth.org/>

614

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

Figure 1. Polyester anterior nares swabs are both comfortable and easy to use. (A,B)

Study participants' (n=35) self reported **(A)** discomfort and **(B)** confidence during self administration of an anterior nares swab at home. **(C)** Boxplot depicting the RT-qPCR crossing threshold values for RNaseP from self-administered anterior nares swabs (ANS) and mid-turbinate (MT) swabs.

Figure 2. Extraction free RT-qPCR set-up and test performance. (A) Assay layout of the EF-RT-qPCR test. One sample is assayed in four wells on a 384 well-plate. Each sample is tested for two probes, in duplicate. RNase P is assayed in each well. **(B)** Mean Ct values for 67 specimens processed by EF-RT-qPCR and extraction-based RT-qPCR. Reactions with no amplification by one preparation protocol are demarcated with red or blue points as indicated.

Figure 3. Addition of Proteinase K improves test performance. (A) Observed percentage of test failures with and without the addition of Proteinase K. **(B)** Archived samples with Ct > 28 reprocessed with either KingFisher Flex Extraction (left), Extraction-Free RT-PCR (middle), or SwabExpress (Extraction-Free RT-PCR + ProK) (right). Colors signify the number of samples and their classifications. **(C)** Box and whisker plots depicting the average delta Ct between replicate wells for SARS-Cov2 positive specimens. Red points indicate outliers. ΔCt values were more consistent upon addition of Proteinase K. **(D)** Mean Ct values of matched specimens run through the automated KingFisher extraction system (left) or using SwabExpress (SE). Specimens which were detected in only one of the two protocols are displayed as green points.

Figure 4. SwabExpress workflow. (1) Anterior nares swabs are collected and (2) transported dry to the lab. Upon receipt, (3) each swab is then hydrated with low-TE buffer, aliquoted into a 96 well plate, and (4) proteinase K is added to every well. (5) The eluted specimens are digested

1
2
3 640 and heat-inactivated in a laboratory oven before (6) they are loaded as the template in a RT-
4
5 641 qPCR reaction. The cost listed includes reagents and consumables.
6
7 642
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60




Figure 1. Polyester anterior nares swabs are both comfortable and easy to use. (A,B) Study participants' (n=35) self reported (A) discomfort and (B) confidence during self administration of an anterior nares swab at home. (C) Boxplot depicting the RT-qPCR crossing threshold values for RNaseP from self-administered anterior nares swabs (ANS) and mid-turbinate (MT) swabs.

147x52mm (300 x 300 DPI)




Figure 2. Extraction free RT-qPCR set-up and test performance. (A) Assay layout of the EF-RT-qPCR test. One sample is assayed in four wells on a 384 well-plate. Each sample is tested for two probes, in duplicate. RNase P is assayed in each well. (B) Mean Ct values for 67 specimens processed by EF-RT-qPCR and extraction-based RT-qPCR. Reactions with no amplification by one preparation protocol are demarcated with red or blue points as indicated.

216x82mm (300 x 300 DPI)




Figure 3. Addition of Proteinase K improves test performance. (A) Observed percentage of test failures with and without the addition of Proteinase K. (B) Archived samples with Ct > 28 reprocessed with either KingFisher Flex Extraction (left), Extraction-Free RT-PCR (middle), or SwabExpress (Extraction-Free RT-PCR + ProK) (right). Colors signify the number of samples and their classifications. (C) Box and whisker plots depicting the average delta Ct between replicate wells for SARS-Cov2 positive specimens. Red points indicate outliers. Δ Ct values were more consistent upon addition of Proteinase K. (D) Mean Ct values of matched specimens run through the automated KingFisher extraction system (left) or using SwabExpress (SE). Specimens which were detected in only one of the two protocols are displayed as green points.

158x129mm (300 x 300 DPI)

SwabExpress: \$5.98 per test




Figure 4. SwabExpress workflow. (1) Anterior nares swabs are collected and (2) transported dry to the lab. Upon receipt, (3) each swab is then hydrated with low-TE buffer, aliquoted into a 96 well plate, and (4) proteinase K is added to every well. (5) The eluted specimens are digested and heat-inactivated in a laboratory oven before (6) they are loaded as the template in a RT-qPCR reaction. The cost listed includes reagents and consumables.

201x91mm (300 x 300 DPI)